

Alioth Disposable Connector Kit

Disposable Design

Strong Compatibility

Stable Performance



⊕ Product Introduction

Alioth disposable straight-through connector and Blind-plug connector sets, which are made from polypropylene, can enable the correct installation of depth filter and meet the requirements of depth filtration process when used with Alidep depth filter and holders.

⊕ Key Features and Benefits

- Good chemical compatibility, with a tolerance range of pH 1-14
- Good pressure resistance and sealing performance
- Single use recommended
- Stable performance, sturdy structure and durable service
- Available in straight-through and blind-plug connectors

⊕ Ordering Information

Order No.	Product Description	Package Available for Ordering
KITORING	3 straight-through connectors,	1=1/PK
	3 blind-plug connectors	
KITORING 6	6 straight-through connectors	1=1/PK

